

KW 14 Suspension tray with lid



When using chemically aggressive liquids (e.g. TICKOPUR J 80 U), removing chemically aggressive or abrasive contaminants or sonicating sample liquids, insert cups are used (indirect sonication). Because nothing may be placed on the bottom of the ultrasonic bath, BANDELIN offers suitable insert cups made of plastic with a lid. Please note that in order to transmit the ultrasound to the insert beaker, the surface tension of the water must be reduced and therefore a cleaning solution (e.g. TICKOPUR R 33 3%) must be used as coupling liquid / contact liquid.

Plastic insert vessel with lid

for stainless steel aggressive chemicals

- for RK 514 / RK 514 H / ZE 514
- inside: L x W x D = 280 x 215 x 145 mm
- temperature resistant: Water up to 80 °C / Acids up to 60 °C
- Polyethylene, unperforated

item number

WL18613

model

KW 14

Products for the electronic industry



manufacturer	BANDELIN
manufacturer item number	613
length	280 mm
width	215 mm
height	145 mm
order unit	1 piece
content unit	1 piece