


DE 100 Perforated lid for ultrasonic bath



When using chemically aggressive liquids (e.g. TICKOPUR J 80 U), removing chemically aggressive or abrasive contaminants or sonicating sample liquids, insert beakers are used (indirect sonication). Because nothing may be placed on the bottom of the ultrasonic bath, BANDELIN offers suitable perforated lids/insert strips to hold the insert beakers. Further advantage: Several insert cups can be sonicated with different liquids at the same time (e.g. removal of dental cements with STAMMOPUR Z in one and cleaning with simultaneous chemical disinfection of dental instruments with STAMMOPUR DR 8 in a second insert cup).

- Stainless steel
- for the insertion of 2 insertion cups
- Suitable for DL 102 H

Item no.

WL20576

Model	DE100
Manufacturer	BANDELIN
Manufacturer article no.	3017
GPSR manufacturer data	BANDELIN electronic GmbH & Co. KG Heinrichstrasse 3-4 DE-12207 Berlin www.bandelin.com
Length	264 mm
Width	132 mm
Height	11 mm
Sales unit	1 piece
Content unit	1 piece
Incl. battery	no