

## DE 510 Perforated lid for ultrasonic bath



When using chemically aggressive liquids (e.g. TICKOPUR J 80 U), removing chemically aggressive or abrasive contaminants or sonicating sample liquids, insert beakers are used (indirect sonication). Because nothing may be placed on the bottom of the ultrasonic bath, BANDELIN offers suitable perforated lids/insert strips to hold the insert beakers. Further advantage: Several insert cups can be sonicated with different liquids at the same time (e.g. removal of dental cements with STAMMOPUR Z in one and cleaning with simultaneous chemical disinfection of dental instruments with STAMMOPUR DR 8 in a second insert cup).

- Stainless steel
- for hanging 4 insert cups.
- Suitable for RK 510/H, RK 512 H, DL 512 H

item number	WL21940
model	DE510
manufacturer	BANDELIN
manufacturer item number	3038
length	329 mm
width	231 mm
height	11 mm
order unit	1 piece
content unit	1 piece