


## DE 514 Perforated lid for ultrasonic bath

---



---

When using chemically aggressive liquids (e.g. TICKOPUR J 80 U), removing chemically aggressive or abrasive contaminants or sonicating sample liquids, insert cups are used (indirect sonication). Because nothing may be placed on the bottom of the ultrasonic bath, BANDELIN offers suitable perforated lids/insert strips to hold the insert beakers. Further advantage: Several insert cups can be sonicated with different liquids at the same time (e.g. removal of dental cements with STAMMOPUR Z in one and cleaning with simultaneous chemical disinfection of dental instruments with STAMMOPUR DR 8 in a second insert cup).

- Stainless steel
- for hanging 4 insert cups.
- Suitable for RK 514/H/BH, DL 514 BH

Item no.

WL21941

Model	DE514
Manufacturer	BANDELIN
Manufacturer article no.	3039
GPSR manufacturer data	BANDELIN electronic GmbH & Co. KG Heinrichstrasse 3-4 DE-12207 Berlin www.bandelin.com
Length	356 mm
Width	291 mm
Height	11 mm
Sales unit	1 piece
Content unit	1 piece