

KW 10-0 Suspension tray



Insert tubs are used for the application of chemically aggressive liquids (e.g. TICKOPUR J 80 U), the removal of chemically aggressive or abrasive contaminants or the sonication of sample liquids (indirect sonication). Because nothing may be placed on the bottom of the ultrasonic bath, BANDELIN offers suitable insert cups made of plastic with a lid. Please note that in order to transmit the ultrasound to the insert beaker, the surface tension of the water must be reduced and therefore a cleaning solution (e.g. TICKOPUR R 33 3%) must be used as coupling liquid / contact liquid.

- Plastic suspension tank with lid
- for stainless steel aggressive chemicals.
- for RK 510 / RK 510 H
- inside: L x W x D = 242 x 182 x 136 mm
- temperature resistant: Water up to 80 °C / Acids up to 60 °C
- Polyethylene, unperforated

Item no.	WL22336
Model	KW 10-0
Manufacturer	BANDELIN
Manufacturer article no.	3053
GPSR manufacturer data	BANDELIN electronic GmbH & Co. KG Heinrichstrasse 3-4 DE-12207 Berlin www.bandelin.com
Length	242 mm
Width	182 mm
Height	135 mm
Sales unit	1 piece
Content unit	1 piece
Incl. battery	no