

SD 04 Insert beaker cpl. for ultrasonic bath



Glass, 400 ml, for insertion into perforated cover/insertion strip with ring GR 04. Insertion beakers are used for the application of chemically aggressive liquids (e.g. TICKOPUR J 80 U), the removal of chemically aggressive or abrasive contaminants or the sonication of sample liquids (indirect sonication). Because nothing may be placed on the bottom of the ultrasonic bath, BANDELIN offers suitable perforated lids/insert strips to hold the insert beakers. Further advantage: Several insert cups can be sonicated with different liquids at the same time (e.g. removal of dental cements with STAMMOPUR Z in one and cleaning with simultaneous chemical disinfection of dental instruments with STAMMOPUR DR 8 in a second insert cup). Please note that in order to transfer the ultrasound to the insert cup, the surface tension of the water must be reduced and therefore a cleaning solution (e.g. TICKOPUR R 33 3%) must be used as coupling liquid / contact liquid.

□ to hang in perforated cover DE08

□ suitable for RK31/H

□ D = 76 mm, H = 110 mm

□ capacity = 600 ml

□ Glass, complete with retaining ring

item number	WL27426
model	SD 04
manufacturer	BANDELIN
manufacturer item number	168
length	75 mm
width	80 mm
height	110 mm
order unit	1 piece
content unit	1 piece

