

## PNG reference app



With PNG Reference, a software app for TAGARNO digital microscopes, you can superimpose self-created inspection images and reference images on the live microscope image. We suggest creating inspection images to verify sample sizes or reference images to sort batches and find defects in products. Once selected, the image is overlaid on the live microscope image to assist you in quality control. With this TAGARNO digital microscope application, you can save your photo along with the microscope settings. When you recall the preset, the microscope automatically changes to the magnification and light settings defined in the preset. This is a much faster and more accurate process than if you had to do it manually. A mouse is required to use this application.

item number	WL47699
model	306010
manufacturer	TAGARNO
manufacturer item number	306010
order unit	1 piece
content unit	1 piece

Products for the electronic industry



ESD safe

no

