

Eye wash station Medidrop II 1x500ml sodium chloridl. 1x250ml BioPhos74



With an eye injury, every second counts. In order to be able to provide immediate first aid, it is highly recommended that you store the eye rinsing bottles in a special wall holder. In addition, make sure that the eyewash bottles are available in a fixed and easily accessible place. The Eye Care Medidrop II eye wash station contains two different eye wash bottles for different applications:

 1×500 ml sodium chloride solution, .9%: for rinsing out foreign bodies such as dirt, dust, chips, wood and metal splinters 1×250 ml BioPhos74 with phosphate-buffered rinsing solution (4.9%) for rinsing out chemical, acid and alkali burns

The eye wash stations from Franz Mensch are dust-tight and space-saving. You will always receive a complete set for wall mounting, and an integrated eye emergency plan and explanatory pictograms support the application.

ready to use and sterile eye wash bottles with ergonomic rinsing head incl. mirror and pictogram for easy handling dust-tight and hygienically space-saving











27.04.2024



Click here for the eye emergency plan!

- Plastic, sterile
- \bullet Content: 1 x 500 ml sodium chloride solution (.9%) and 1 x 250 ml phosphate-buffered rinsing solution (4.9%)
- ready to hand, ready to use, dust-tight and space-saving
- easy mounting and handling
- integrated eye emergency plan

item number	WL50190
model	50016
manufacturer	FRANZ MENSCH
manufacturer item number	50016
order unit	1 piece
content unit	1 piece







