

## Tickopur Neutral Cleaner Concentrate, R30 / 2 litres



### Neutral cleaner

For particularly gentle cleaning and degreasing of dismantled and undismantled machines, tool, glass and ceramic parts in industry, trade and laboratories.

TICKOPUR R 30 is specially designed for cleaning in SONOREX ultrasonic units, but can also be used for immersion, wiping and high-pressure cleaning. When used properly, a fast and gentle cleaning is guaranteed.

### TICKOPUR R 30 removes

light grinding, polishing and lapping residues, oil- and grease-like soiling, soot etc. from parts and equipment made of metal, also light metal, glass, ceramics, plastic, rubber etc.

- Application concentration
- Application with ultrasound
- Dosage: 1 - 5
- sonication time: 1 - 10 minutes

- Temperature: 20 - 80 °C

item number	WL18023
model	R30 / 2 Liter
manufacturer	Dr. H. Stamm
manufacturer item number	879
order unit	1 piece
content unit	1 piece