

## HygoClean Alcohol-based rapid disinfection, 10 litres, for LMP / Curacid Aktiv Plus



Important information: This product is classified as a hazardous material. We transport dangerous goods only within the German-speaking countries and only to addresses accessible by rail or road. Please also refer to our General Terms and Conditions (Section V.1).

Rapid disinfection for surfaces: Especially for food areas This highly effective surface disinfection has been specially developed for use in food areas. Because clean surfaces and working materials are particularly important wherever food is handled and food is prepared. This is a ready-to-use disinfection solution that is ready for use immediately after opening. It is purely alcoholic and free of anything that could contaminate food - without colorants or fragrances. The rapid disinfection is ideal for all surfaces, machines, knives or other work equipment in the food industry and in kitchens. In contrast to cleaning, where germs are only washed off, disinfection kills the germs. In hygiene-sensitive areas, disinfection therefore provides a more reliable germ stop than cleaning - for optimum,

hygienic results. The disinfectant solution can simply be used directly with a cloth for wiping. Alternatively, the solution can also be sprayed onto the surface as spray disinfection with the appropriate attachment.

extremely effective against bacteria and fungi food-safe fast-acting, after only one minute effective residue-free evaporates very quickly and dries completely degreases and disinfected for wiping: Simply pour the finished solution onto a cloth and wipe the dry surface and rinse with water (not suitable for acrylic/plexiglass) as a spray: suitable attachment available, simply spray the solution onto the surface with a cloth and rinse with water Effective spectrum Rapid disinfectant: Bactericidal and levurocidal effect according to EN 1276, EN 165, EN 13967. Valid for Europe-wide use in the food sector. Meets the requirements of the European food law (HACCP). Two sizes available: large 1 litre canister for high demand or practical 1 litre bottle Spray attachment for 1 litre bottle available Canister tap for 1 litre canister for easy removal available

The quick disinfection for surfaces can be used in combination with the HygoWipe system (consisting of cloth and bucket). Thus you have an extremely effective and practical disinfection system, with which you also work very economically and economically.

- ready-to-use alcoholic solution
- immediately ready for use - simply apply to a cloth (e.g. Hygowipe), wipe dry surface and rinse with water (not for acrylic/ plexiglass), without fragrances or dyes
- evaporates very quickly and dries completely, degreases and disinfects
- Use biocidal products carefully. Always read label and product information before use

<b>Item no.</b>	<b>WL49698</b>
Model	HygoClean-31565
Manufacturer	FRANZ MENSCH
Manufacturer article no.	31565
GPSR manufacturer data	Franz Mensch GmbH Werner-von-Siemens-Straße 2 DE-86807 Buchloe www.franz-mensch.de
Sales unit	1 jerrycan
Content unit	10 l
ESD compliant	no